



SOP GREEN KLEAN®
Chlorinated Disinfecting Tablets



DIRECTIONS FOR USE GREEN KLEAN CHLORINATED DISINFECTING TABLETS

Container Size	Type Solution	Number of Tablets
5 Gallon Buddy Jug	EPA Sanitizer	1
Mop Bucket – Up to 5 Gallons	EPA Sanitizer	1
Third Sink – 5 Gallons	EPA Sanitizer	1
Quart Spray Bottle	EPA Disinfectant Kills Norovirus in 1 minute C Diff in 10 minutes	1
Quart Spray Bottle	EPA Disinfectant Kills C Diff in 4 minutes	2
Quart Spray Bottle	Grout Cleaner / Mold and Mildew Stain Remover	2
Electrostatic Sprayer Quart Tank	EPA Disinfectant @ 2153 ppm Effective against 31 pathogens	1

*SOLUTION RETAINS EFFICACY FOR 3 DAYS WHEN ADDED TO WATER AND KEPT IN A CLOSED CONTAINER.



Distributed By:



VICTORY
INNOVATIONS C9

VP200ESK

BETTER VISIBILITY:
LED HEAD LIGHT

INCREASED DURABILITY:
GLASS-FILLED HOUSING



ERGONOMIC
HANDLE WITH LOCK
FOR COMFORT AND
SAFETY

LONG LASTING BATTERY VP20
(4 hr. Run Time)

VERSATILITY:
ADJUSTABLE 3-IN-1
NOZZLE

DOUBLE CHARGE
TECHNOLOGY FOR
MAXIMUM COVERAGE

1 LITER EASY
FILL TANK

HEAVY DUTY
KIT BOX



PROFESSIONAL CORDLESS ELECTROSTATIC HANDHELD SPRAYER

PRODUCT SPECIFICATIONS

Model #	VP200ESK	Weight (No Chemical)	3.8 lbs.
Power	Cordless	Weight (Full)	5.9 lbs.
Tank Size	33.8 ounces / 1 L	Optimum Spray Range	2-3 ft.

	NOZZLE SETTING 1	NOZZLE SETTING 2	NOZZLE SETTING 3
Nozzle Shape	Full Cone	Full Cone	120° Fan
Particle Size	40 microns	80 microns	110 microns
Flow Rate	3.1 ounce/min	3.8 ounce/min	10.5 ounce/min
Run Time per Tank	11.3 min.	9.3 min	3.33 min
Coverage per Tank*	2,800 sq ft	2,550 sq ft	1,700 sq ft
Tanks per Charge	21.2	25.8	72.1

*Per internal testing determined by flow rate and particle size

ACCESSORIES:



VP72 / VP74
Extension Wand -
12" or 24"



VP30 33.8oz
Tank with Cap



VP49 Nozzle
Wrench



VP50 3-in-1
Nozzle



VP20B 16.8V
6800mAh Battery
(8 hr. run time)



VP91 Carry
Strap

FC CE RoHS NOM



VICTORY
INNOVATIONS CO

VP300ES

PORTABILITY:
— EASY CARRY HANDLE

2.25 GALLON EASY FILL
TANK



VERSATILITY:
ADJUSTABLE 3-IN-1
NOZZLE

COMFORTABLE
CARRYING STRAPS

DOUBLE CHARGE
TECHNOLOGY FOR
MAXIMUM COVERAGE

INCREASED DURABILITY:
GLASS-FILLED HOUSING

ERGONOMIC
HANDLE WITH LOCK
FOR COMFORT AND
SAFETY

QUICK RELEASE
VALVE

LONG LASTING BATTERY VP20A
(4 hr. Run Time)

— 4 FT HOSE

PROFESSIONAL CORDLESS ELECTROSTATIC BACKPACK SPRAYER

PRODUCT SPECIFICATIONS

Model #	VP300ES	Weight (No Chemical)	10 lbs.
Power	Cordless	Weight (Full)	28.8 lbs.
Tank Size	2.25 gallons / 8.52 L	Optimum Spray Range	4 - 6 ft.

	NOZZLE SETTING 1	NOZZLE SETTING 2	NOZZLE SETTING 3
Nozzle Shape	Full Cone	Full Cone	120° Fan
Particle Size	40 microns	80 microns	110 microns
Flow Rate	3.4 ounce/min	5.1 ounce/min	13.4 ounce/min
Run Time per Tank	84 min	56.5 min	21.8 min
Coverage per Tank*	23,000 sq ft	20,600 sq ft	13,100 sq ft
Tanks per Charge	2.9	4.2	11

*Per internal testing determined by flow rate and particle size

ACCESSORIES:



VP72 / VP74
Extension Wand -
12" or 24"



VP31 2.25 Gallon
Tank with Cap



VP49 Nozzle
Wrench



VP50 3-in-1
Nozzle



VP20B 16.8V
6800mAh Battery
(8hr. run time)



VP10 16.8V
Battery Charger

FC ETL CE RoHS NOM



DISINFECT ALL COMMERCIAL FACILITIES FROM ENTRY TO EXIT

EPA REGISTERED USES: For use in amusement parks, breweries, beverage and food processing plants, food processing and canning plants, bottling plants, fish processing plants, meat and poultry processing plants, schools, universities, hospitals, neo-natal units, physicians/pediatricians offices, office buildings, nursing homes, elder care centers, child care centers, daycares, shelters, cafeterias, diners, eating establishments, institutional dining establishments, hotel restaurants, food service, stores, shops, camp sites, casinos, movie theaters, industrial facilities, military installations, licensed care facilities, institutions, catering, kitchens, vending machines, Intensive Care Unit, ICU operating rooms, emergency waiting rooms, dental facilities, gyms, locker room facilities, health clubs, restrooms, Tattoo parlors, manicure pedicure salons, commercial laundries, and prisons. This product is for use in hospitals, medical and dental offices and clinics, operating rooms, isolation wards, and medical research facilities. Cosmetic, pharmaceutical, medical device manufacturing facilities, biotechnology firms, pharmacies, and compounding pharmacies.



6.55 Gram Tablet

PATHOGEN	MINIMUM DOSE REQUIRED	MINIMUM CONTACT TIME REQUIRED
Food Contact Sanitizer Claims		
Staphylococcus aureus	100 ppm	1 minute
Salmonella enterica	100 ppm	1 minute
Disinfection Claims – bacteria		
Staphylococcus aureus	a) 1076 ppm b) 4306 ppm	10 minutes 4 minutes
Staphylococcus aureus – methicillin resistant (MRSA) & glycopeptide-resistant (GRSA)	a) 1076 ppm b) 4306 ppm	10 minutes 4 minutes
Staphylococcus epidermidis	1076 ppm	10 minutes
Salmonella enterica	a) 1076 ppm b) 4306 ppm	10 minutes 4 minutes
Pseudomonas aeruginosa	a) 1076 ppm b) 4306 ppm	10 minutes 4 minutes
	4306 ppm	4 minutes
Streptococcus pneumoniae		
Escherichia coli O157-H7	1076 ppm	10 minutes
Acinetobacter baumannii	4306 ppm	4 minutes
Vancomycin resistant Enterococcus faecalis	a) 1076 ppm b) 4306 ppm	10 minutes 4 minutes
Carbapenem resistant Klebsiella pneumoniae	4306 ppm	4 minutes
Klebsiella pneumoniae	1076 ppm	10 minutes

Effective against gram positive and gram negative bacteria.

PATHOGEN	MINIMUM DOSE REQUIRED	MINIMUM CONTACT TIME REQUIRED
Virucidal Claims		
Respiratory syncytial virus	538 ppm	10 minutes
Rhinovirus Type 14	1076 ppm	10 minutes
Influenza H1N1	1076 ppm	10 minutes
Human Immunodeficiency Virus Type 1 (HIV-1)	a) 1076 ppm b) 4306 ppm	10 minutes 1 minute
Hepatitis A and B virus	a) 1076 ppm b) 4306 ppm	10 minutes 1 minute
Hepatitis C virus	4306 ppm	1 minute
Coxsackievirus	4306 ppm	1 minute
Influenza	4306 ppm	1 minute
Norovirus	2153 ppm	1 minute
Poliovirus Type 1	1076 ppm	10 minutes
Fungicidal/Yeastcidal Claims		
Aspergillus fumigatus	4306 ppm	1 minute
Candida albicans	4306 ppm	4 minutes
Trichophyton mentogrophytes	1076 ppm	10 minutes
Herpes simplex virus type 1	1076 ppm	10 minutes
Clostridium difficile Claims		
Clostridium difficile spores	a) 2153 ppm b) 4306 ppm	10 minutes 4 minutes
Mycobactericidal Claims		
Mycobacterium bovis (TB)	5382 ppm	4 minutes

SANITIZING DIRECTIONS

TABLET SIZE 6.55 grams	TABLETS PER 5 GALLONS OF WATER
100 PPM	1

DISINFECTION MIXING DIRECTIONS

TABLET SIZE 6.55 grams	TABLETS PER QUART OF WATER	TABLETS PER GALLON OF WATER
538 PPM	0	1
1076 PPM	0	2
2153 PPM	1	4
4306 PPM	2	8
5382 PPM	0	10



See page 2 for Directions for use

SCHOOLS DEMAND SUCCESS IN CONTROLLING

☑ **ODORS** ☑ **MOLD**

&

☑ **ICE MACHINE**
Pathogens

They are potential hazards to plan for, deal with, and overcome—proactively.....



INSPIRED TEC HAS YOUR SOLUTION

Our air and surface purification technology is a **NON-CHEMICAL—ACTIVE strategy** for eliminating **ODORS, MOLD,** and effectively purifying your



A Multi Use Machine...

Inspired TEC's PCO (Photo-catalytic Oxidation) technology, using an enhanced ionization output, represents the newest generation of proactive air purification. Our unit is designed for any application, and protects indoor environments with 24/7 active air and surface purification—PLUS.....

Ice Machines achieve optimal results. **CLEAR** technology has been tested and proven to significantly reduce Listeria, biofilms, E-Coli, Aspergillus, Salmonella, Norwalk, Mold and other pathogens in Ice Machines and coolers...PLUS.

